

Instructions for use

Eurogloves Disposable Gloves, Vinyl, (Soft) Nitrile & Latex.



All Eurogloves disposable gloves are in conformity with EU Community Legislation Medical Devices Regulation (EU) 2017/745 and Personal Protective Equipment (PPE) Regulation (EU) 2016/425.

The gloves are distributed within a dispenser. In case the dispenser has been dismantled and the gloves are sold separately, the distributor must ensure that each unit is accompanied with the instructions for use.

All Eurogloves are non-sterile and designed for single-use only. The gloves are intended to be worn on the hand during medical examinations and non-invasive procedures to provide a barrier between the patient and the medical professional.

The gloves meet the applicable EN/EN-ISO standards displayed on each dispenser.

All Eurogloves disposable gloves are classified as:

- Personal Protective Equipment (PPE) Category III, according to PPE Regulation (EU) 2016/425.
- Medical Device Class I, according to Medical Device Regulation (EU) 2017/745 Annex VIII per rules 1 and 5.

Icon Explanation



Single-use only. Not designed to be reused.



Product is safe for food contact.



Accepted Quality Limit: Percentage of products that could deviate from standards.



Carton dispenser can be recycled.



Store product out of rain.



Product is classified as a medical device.



Store product between 5°C and 30°C.



Store product away from sunlight.



Tested for resistance to penetration by bloodborne pathogens according to ASTM.



Number of notified body responsible for ongoing EU type examination.



EN-ISO 374-1:2016 +A1:2018 level of chemical penetration resistance based on testing according to EN 16523-1:2015 against the specified test chemicals.

ABCDEFGHIJKLMNOST

EN-ISO 374-1: 2016 +A1:2018 Test Chemicals.

- | | |
|-----------------------|-----------------------------|
| A= Methanol | J = n-Heptane |
| B= Acetone | K = Sodium hydroxide 40 % |
| C= Acetonitrile | L = Sulphuric acid 96 % |
| D= Dichloromethane | M = Nitric acid 65 % |
| E = Carbon disulphide | N = Acetic acid 99 % |
| F = Toluene | O = Ammonium hydroxide 25 % |
| G = Diethylamine | P = Hydrogen peroxide 30 % |
| H = Tetrahydrofuran | S= Hydrofluoric acid 40 % |
| I = Ethyl acetate | T= Formaldehyde 37% |

Standards Explanation

EN-ISO 374-1: 2016 +A1:2018	Specifies performance requirements against chemical penetration.
Levels of chemical penetration resistance: <i>Type A= Chemical breakthrough time >30 minutes against at least 6 test chemicals.</i> <i>Type B= Chemical breakthrough time >30 minutes against at least 6 test chemicals.</i> <i>Type C= Chemical breakthrough time >30 minutes against at least 1 test chemical.</i>	
EN-ISO 374-2: 2019	Specifies a test method for the penetration resistance of chemicals and micro-organisms.
EN-ISO 374-4: 2019	Specifies a test method for the penetration resistance of chemicals with continuous contact.
EN-ISO 374-5: 2016	Specifies performance requirements and test methods against micro-organisms.
EN-455-1:2020	Specifies performance requirements and test methods for freedom from holes.
EN-455-2:2015	Specifies performance requirements and test methods for physical properties.
EN-455-3:2015	Specifies performance requirements for the evaluation of biological safety.
EN-455-4:2009	Specifies performance requirements for shelf life.
EN-ISO 21420:2020	Specifies general requirements and relevant test procedures.

Precautions for use

Before use, inspect the gloves for any defects, e.g. holes and tears. Do not use defective gloves.
 If in doubt, do not use the gloves, get a new pair.

Storage Instructions

Store this product in a dry and dust free environment between 5°C and 30 °C. Keep product in a dark place, away from (ultraviolet) light sources.

Hazardous Ingredients

Some variants of Eurogloves contain ingredients which are known to possibly cause skin irritation or allergic reactions. Always check warning information on dispenser before use.

Disposal Instructions

Used gloves may be contaminated with infectious or other hazardous materials. Dispose of according to local authority regulations.