



Traverse Rescue Basket Stretchers



When you're called upon to perform remote or technical rescue, you need to have equipment on hand that you can trust.

FERNO offers a wide range of reliable and robust **Rescue Basket Stretchers** from Traverse Rescue (a FERNO Group Company) that are easy to handle, safe for patient and rescuer, and go wherever rescue is required.

From premium options like the **Titan Series** to the capable-yet-cost-effective **Gazelle Series**, FERNO has a basket stretcher that will meet your needs.

- **Reliable and durable**
basket stretchers for all budgets
- **Wide variety of options**
for the versatility to tackle any situation
- **Incredible strength and durability** to hold up to the demands of rescue
- **Designed for safety and comfort** for patients and operators



Split-apart models facilitate carrying and storage and are ideal for helicopter rescue



Tapered models allow for easier access into confined areas



StratLoad™ attachment points, featured as standard on several models, enable fast and safe attachment of litter bridle

Specifications



	Advantage	Bradco	Gazelle (regular)	Spartan (titanium)	Titan (regular)	TRS
LENGTH	85 in 216 cm	80 in 203 cm	82.5 in 210 cm	86 in 218 cm	83.5 in 212 cm	79.5 in 202 cm
WIDTH	24 in 60 cm	25.2 in 64 cm	24.5 in 62 cm	24.5 in 58 cm	23 in 58 cm	33 in 84 cm
HEIGHT	7.5 in 19 cm	8.5 in 21.6 cm	7 in 17 cm	7.5 in 18.41 cm	7.25 in 18 cm	N/A
WEIGHT	33 lb 15 kg	46 lb 20.8 kg	33 lb 15 kg	22 lb 10 kg	30 lb 14 kg	18 lb 8 kg
LOAD LIMIT	2500 lb 11 kN	2500 lb 11 kN	900 lb 4 kN	2500 lb 11 kN	2500 lb 11.3 kN	550 lb 250 kg (vert) 771 lb 350 kg (horiz)

Features

Advantage

- Fully-internal stainless-steel frame provides 2,500-lb strength rating and a smooth, single-plane sliding surface
- Durable HDPE shell is attached to the top tube, providing numerous options for securing patients to the exposed secondary rail
- Cut-outs in shell provide multiple sections of exposed top rail for handling
- Closed-cell foam mattress is impervious to fluids
- StratLoad attachment points featured as standard

Bradco

- Well-suited for any type of rescue, including industrial applications
- Designed for horizontal and vertical use
- Heavy-duty, all-aluminum frame protects patient and has a 2,500-lb strength rating
- Spring-suspended patient support and restraint system increases patient safety and comfort
- Double-flap cover allows patient to be transported securely in almost any position
- Polyethylene shell is resistant to chemicals and glides smoothly over various types of terrain

Gazelle Series

- Cost-effective and capable basket stretchers tested to and beyond a 900-lb load capacity
- Incredibly versatile and available in Regular, Tapered, and Confined Space models
- All-carbon-steel MIG-welded frame for exceptional durability and strength
- Polyethylene netting is comfortable and allows water to pass through easily

Spartan

- Strong stainless-steel frame and hard polyethylene shell offer excellent durability
- Split-Apart design facilitates carrying and storage
- Fully-exposed top rail offers easy handling
- Recessed sliding surface is capable of handling rough terrain
- StratLoad attachment points featured as standard
- Available in stainless steel and titanium

Titan Series

- Designed to perform in the most technical rescue environments
- Offers 2,500-lb load capacity
- Available in stainless steel and titanium
- Incredibly versatile and offered in One-Piece or Split, Regular or Tapered models
- Polyethylene netting is comfortable and allows water to pass through easily
- StratLoad attachment points featured as standard

TRS

- Designed for horizontal and vertical evacuations
- Wrap-around design offers comfort and security to patients
- High-density polyethylene construction is durable and protects patient
- Built in restraint system is secure and quick to apply

Accessories



Soft Eye Litter Bridle
No. 0313898



Model 418-1 Adjustable Lifting Bridle
No. 0312390



Model 132 Manta Rescue Aid
No. 0382888 Orange
No. 0382969 Burgundy
No. 0383034 Olive Drab Green



Hard Eye Litter Bridle
No. 0313897



Model 471 Flotation Collar
No. 0313943

NOT PICTURED

Carrying/Storage Cases

Contact your sales representative for more information