



NOBAMED Paul Danz AG

KOMPRESSEN-steril

REF 712105

Product Description and Purpose

Gauze swabs, **12-ply**, made of 17 thread dressing gauze, **size 5 cm x 5 cm**, individually **packed per 1 piece** in sterile packaging, raw edges folded back. The gauze swabs are used in primary wound care, for wound cleansing, wound covering, and for the absorption of blood and wound fluid. They are for single use only and have been packaged in a way that their folded back points upward, thus allowing an easy removal.

Contra-Indications

No contra-indications known. The product should not be used in the case of a known allergy to the ingredient listed below.

Composition

Cotton

Normative and Legal Requirements

KOMPRESSEN-steril (sterile swabs) are medical products according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and are classified as class Is, rule 4.

The gauze swabs comply with the requirements of the DIN EN 14079.

Sterilization of the product complies with DIN EN ISO 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	paper-paper-packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols

