

ri-gital®
Easy to use, robust and precise.

ri-gital® Digital Thermometer

Clinical thermometer in the well-known Riester quality. Can be used for oral, rectal and under-arm measurement of body temperature. Digital measurement is particularly accurate thanks to the highly sensitive temperature sensor.

- › 3-digit display, memory and signal tone
- › Water-proof
- › Measurement range: 32 – 43.9°C
- › Accuracy: +/- 0.1°
- › Automatic switch-off after 10 min
- › Includes long-life button cell batteries
- › Comes in plastic storage box
- › Available with disposable probe covers (not included)

No. 1851	ri-gital® digital thermometer with button cell battery
No. 12703	Disposable probe covers ri-gital® , pack of 40 pcs.
No. 12706	Disposable probe covers ri-gital® , pack of 1000 pcs.

Please visit **riester.de** for more information