



**NOBAMED Paul Danz AG**

## LIEGENPAPIER

**REF 650749**

### **Product Description and Purpose**

The medical paper roll (49cm x 150m) has a surface weight of 35 g/m<sup>2</sup> and has a perforation after each 36 cm. It is used for hygienic covering of the examination couch.

### **Composition**

Cellulose

### **Normative and Legal Requirements**

Liegenpapier (medical paper rolls) are subject to the Consumer Goods Ordinance.

The product does not contain dangerous toxic substances according to REACH.

It has DIN EN ISO 15223-1 labels on all its packaging.

### **Packaging**

Primary packaging: paper wrap

Secondary packaging: carton made of cellulose

### **Storage**

To be stored in a dry and dust-free environment

### **The product bears following labels and symbols:**

