

ri-thermo[®] tymPRO+

Clinical grade, tympanic ear thermometer with Bluetooth[®] connectivity

- › Clinical grade readings in 1 - 2 seconds
- › Guide indicator to ensure proper probe tip placement
- › Applicable for all patient types
- › Riester, single use, latex and BPA free disposable probe covers.
- › Aid prevention of cross infection with contactless one hand operation
- › Minimum of 3000 measurements with one set of AA batteries
- › Available with Bluetooth for connectivity with Riester **RVS-100** patient monitor, allowing seamless integration with the bedside device
- › Optional calibration kit available
- › CE / FDA approved



Automatic On/Off

Dual-colour, backlit display for low light situations

Switchable from Celsius to Fahrenheit



The Riester **ri-thermo[®] tymPRO+** is a clinical grade tympanic thermometer with guiding indicator and advanced measuring technology. The unique Riester probe cover with dispenser box provides an efficient workflow, minimising cross contamination. Designed for the professional.

Measurement range:	32.0 °C (89.6 °F) – 43.0 °C (109.4 °F)
Accuracy:	±0.2 °C (±0.4 °F) for the range of 35.0 °C to 42.0 °C (95 °F to 107.6 °F) ±0.3 °C (±0.5 °F) for the range of < 35.0 °C (95 °F) or > 42.0 °C (107.6 °F)
Operating Condition:	10 °C to 40 °C (50 °F to 104 °F), 95 % RH or less non-condensing
Storage Condition:	-20 °C to 60 °C (-4 °F to 140 °F), 95% RH or less non-condensing
Resolution:	0.1 °C / 0.1 °F
Measurement Unit:	°C (Default) / °F
Display Function:	Hi / Lo Indicator
Power Source:	1.5 V AA x 2 non rechargeable
Battery Life:	At least 3000 measurements
Memory Capacity:	10 measurements
Power saving:	60 seconds
Communication:	Bluetooth Smart (optional)
Dimension / Weight:	150 (L) x 42 (W) x 56 (H) mm / 80 g
Connectivity:	Optional Bluetooth enabled
No. 1835-BT	Tympanic thermometer ri-thermo[®] tymPRO+ Bluetooth
No. 1835	Tympanic thermometer ri-thermo[®] tymPRO
No. 12634	Probe cover for ri-thermo[®] tymPRO 1 rack / 10 pcs.
No. 12640	Probe cover for ri-thermo[®] tymPRO 30 racks / 300 pcs.
No. 12953	Calibration kit for ri-thermo[®] tymPRO
No. 12911	Anti-theft cable, length 50cm, black, keyring at both end for ri-thermo[®] tymPRO
No. 13347	Bluetooth connectivity: Bluetooth USB-Dongle for RVS-100

Please visit **riester.de** for more information