

ri-clip®

Vein Compressor



Taking blood should be as fast and easy as possible - and this begins with compressing the veins. Discover a new product from **Riester** that makes this task much easier. Because when you put the cuff in place, you have always got one hand free. The highquality plastic (free of all metals) is light as a feather and works like a dream. It makes uncomfortable pinching of the skin or small hairs a thing of the past. And all this at a highly attractive price and in a range of colours. .

- › Simple tension and release clip
- › Aids easy visualisation of the vein for blood draw
- › Hygienic, can be disinfected with conventional disinfectants or autoclaved at 120 °C / 1 bar
- › Available in latex-free, avoiding unpleasant skin reactions
- › Avoids unpleasant catching of small hairs and pinching of skin

No. 5000 ri-clip®

No. LF5000 ri-clip®, latex-free

900521 Rev. E 2023-01-23

Please visit **riester.de** for more information