

AlcoTrue® M

Professional breathalyser for rapid accurate determination of breathalcohol concentration for medical use



» **PRODUCTS** that care
for you «



Easy to use, rapid operation «

AlcoTrue® M – professional breathalyser for rapid accurate determination of breathalcohol concentration in exhaled air. AlcoTrue® M is for diagnostic use, e.g. hospital use, rehabilitation clinics.




Key features

- Identical in construction with the AlcoTrue® P used by police
- Professional **bluepoint®** fuel cell technology
- High precision intelligent microprocessor sampling system
- Automatic self-test after switch on
- Rapid, precise and reliable measurements, even at sub-zero ambient conditions
- Intuitive user interface (self explanatory)
- Ready to test < 6 s, even temperatures below zero
- Results are shown on a colour display (screen)
- Compact, ergonomic design
- Low maintenance and spare part costs
- PC software
- 2 years warranty



bluepoint® MEDICAL sales partner



AlcoTrue® M is certified as class I with measurement function product.

For further information, please visit
www.bluepoint-medical.com

CE 0482

bluepoint MEDICAL